


About Product :-

Product Name	Capping Reagent B 1
Product Combination	40% Acetic Anhydride / Acetonitrile
Catalog#	DN-3306-C
Applications:	DNA/RNA Synthesis
Synonyms	CAP B-1 Capping Reagent B1 (40% Acetic Anhydride in Acetonitrile)
Hazardous Category	
UN Number	2924

Technical Specification

Acetic Anhydride	
Appearance/ Color	Clear and Colorless liquid
Purity by GC	Purity: $\geq 99.9\%$
Boiling Point Range	138°C to 140°C
Acetonitrile	
Appearance/ Color	Clear and Colorless liquid
Purity by GC	Purity: $\geq 99.9\%$
Moisture Content by KF	Moisture Content : <10 ppm
Boiling Point Range	81°C to 82°C
Specification of Cap-B 1 Reagent	
Appearance/ Color	Clear and Colorless liquid
Purity	Acetic Anhydride: 38.0% to 42.0%
Filtered by	0.2 micron GFF
Density	0.9g/cm ³
Storage	Less than 20°C
Shelf Life:	6 months (Under strict storage conditions)

Trade Information

Delivery Time	Immediately
Sample Policy	Paid Sample
Pack Size	Amber Glass Bottle : 2.5 Litre & 4L Litre SS Safety Barrel : 25 Litre & 50 Litre
Packaging Details	UN Approved Packaging
Supply Ability	1000L-1200L per month