


About Product :-

Product Name	20% DEA in Acetonitrile
Product Combination	20% Diethylamine in Anhydrous Acetonitrile
Catalog#	DN-3306-E
Applications:	DNA/RNA Synthesis
Synonyms	20% DEA in Acetonitrile
Hazardous Category	
UN Number	2924

Technical Specification

Diethylamine	
Appearance/ Color	Clear and Colorless liquid
Purity by GC	Purity: $\geq 99.9\%$
Boiling Point Range	54°C to 56°C
Acetonitrile	
Appearance/ Color	Clear and Colorless liquid
Purity by GC	Purity: $\geq 99.9\%$
Moisture Content by KF	Moisture Content : <10 ppm
Boiling Point Range	81°C to 82°C
Specification of Diethylamine Solution	
Appearance/ Color	Clear and Colorless liquid
Purity	Diethylamine: 19.0% to 21.0%
Filtered by	0.2 micron GFF
Storage	Less than 20°C
Shelf Life:	6 months (Under strict storage conditions)

Trade Information

Delivery Time	Immediately
Sample Policy	Paid Sample
Pack Size	Amber Glass Bottle : 2.5 Litre & 4L Litre SS Safety Barrel : 25 Litre & 50 Litre
Packaging Details	UN Approved Packaging
Supply Ability	1000L-1200L per month